



# Donor Donkey Serum

Catalog # 100-151

## PRODUCT DESCRIPTION:

Donor donkey serum is collected aseptically, from normal healthy donor animals.

- Sterile-filtered

## INSTRUCTIONS FOR USE

### Receipt and Storage:

Ships on dry ice. Store frozen at -20 to -10°C

### Thawing Instructions:

When preparing donkey serum for use, it should be allowed to thaw in a refrigerator at a temperature of 2-8°C overnight to minimize thermodynamic stress to the material.

### Appearance:

Straw-like or amber color. It could be slightly cloudy and may contain tiny precipitates.

### Stability:

If stored frozen (-20 to -10°C), product expiration is 5 years from manufacture.

### Preparation and Use:

We recommend that you minimize the number of freeze/thaw cycles the product will undergo by aliquoting the serum into useful volumes at time of storage.

Microbiological Testing	
Bacteria	No Growth
Fungi	No Growth
Physical Testing	
Osmolality	Report
pH	Report
Biochemistry	
Hemoglobin	
Protein, Total	
Storage	Store frozen at ≤ -10°C
Form	Liquid
Expiration	Five years from the date of manufacture

For Cell Culture, further Manufacturing and Research use only. Not for direct Therapeutic use.